
TAPE CORRECTION SOX

TAPE CORRECTION LOOPS

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Effect

- + Combat misalignment and correct malpositioning of toes or foot, depending on the model.
- + Stabilise the toes, mobilise joints, promote a natural rolling movement of the foot, reduce pain and can prevent arthrosis.
- + Suitable for wearing all day and gentle on the skin like normal socks/stockings; discreet can be worn without restrictions and suitable for the workplace, leisure and when exercising
- + To be worn like normal socks/stockings (TAPE SOX) / worn under normal socks or barefoot in summer shoes (TAPE LOOPS)
- + Mobilisation = relief + regeneration

Correction tapes or socks with integrated tapes – they vary according to the model and aim to counteract typical malpositioning.

TAPE SOX and TAPE LOOPS combine the benefits of conservative methods of treatment, such as physiotherapy, with kinesiotaping to counteract malpositioning of the foot or toes. They are free from side-effects. They can be used long term for problems that have functional causes. These include: bunions, flat feet, splayfeet, hammer toes.

Taping Sox and Loops can be used to treat other complaints caused by these misalignments. Along with long-term functional treatment, the position of the foot/toes is immediately altered – for an aesthetic appearance. They are easy to use and suitable for anyone suffering from one of the corresponding misalignments.

Application areas

Preventative: The COMPRESSANA TAPE SOX and TAPE LOOPS can be worn as a precautionary measure if you have an existing or developing misalignment. It is possible and generally unproblematic to combine them with inserts.


Post-operatively: After an operation for bunions or hammer toes, taping socks can help maintain a successful operation result by countering local misalignment.

The TAPE SOX can also be worn at night for gentle, passive correction. The corrective socks offer relief and support for flat feet and splay feet.

TYPE
HALLUX
VALGUS

Bunions (Hallux valgus)

Bunions are the most commonly occurring and treated cause of toe misalignment. Bunions cause the big toe to be malpositioned towards the outer edge of the toe. Women are more prone to this type of foot deformation than men. If left untreated, bunions can cause hammer toes.

TYPE
HALLUX
VALGUS
+PLUS

Hallux valgus +Plus

Advanced hallux valgus can no longer be permanently counteracted by just correcting the big toes. In addition to the toe position, Hallux valgus +Plus additionally corrects the heel position and thereby counteracts the hallux valgus. The combination of the Hallux Valgus Tape Socks and the Pes Valgus Socks makes it possible to counteract further heel cave and the worsening of the hallux valgus tendency. Hallux valgus +Plus is suitable for people with pronounced hallux valgus and splayed feet in combination with pes valgus..

TYPE
HAMMER
TOE

Hammer toes/claw toes

A hammer toe is a toe (usually the second toe) that is stiffened and bent upwards, resembling the shape of a hammer. This malpositioning becomes increasingly painful. Again, women are about five times more likely to suffer hammer toes than men. In the early stages the toes can often be moved back to their normal position. As it progresses, the toes become stiffened and arthrosis develops. Hammer toes and bunions often occur together.

TYPE
FALLEN
ARCH

Flat feet, splay foot (Fallen arch)

If the heel turns inwards and forwards when weight is placed on it, the longitudinal arch becomes flattened. This is termed low arch and flat foot. This causes instability of the outer edge of the heel, making it turn inward. When the body weight is on the foot, there is also pressure on the longitudinal arch. If left untreated, this can cause flat foot. If the shoe sole is worn down towards the inside in the heel and forefoot areas, this indicates overpronation due to fallen arches/ flat foot.

Along with fallen arches and flat feet, TAPE SOX and TAPE LOOPS can be used to treat other foot conditions such as hallux valgus, splay foot, plantar fasciitis and tibialis posterior weakness.

TAPE SOX instructions for use

Bunions / hammer toes: When putting on the socks, please note the separate toe areas. First, gather up the sock and position the big toe in the correct compartment. Do the same for hammer toe corrective socks, putting the second toe in the little box.

The taping band for the hammer toe must be pulled over the middle bone towards the base joint, so that the pull band is over the raised/bent part of the toe (*see image*). The hammer toe tape must not be in front of the raised part as there is a significant risk that this could worsen the misalignment.

Next, position the other toes in the remaining toe boxes and then pull the Taping Sock towards the heel, making sure the pull ties are in the right place. To ensure correct functioning, the taping bands must be under tension with the toes fully in the toe boxes and as taut and crease-free as possible. If the sock is too loose, the corrective function is limited. In this case please choose a smaller size.

Hallux valgus +Plus: Along with a correct fit of the big toes in the toe box, please ensure the correct position of the tape bands in accordance with the following instructions to Fallen arch.

Fallen arch / pronation Control: When putting the socks on, make sure you have the correct socks for right and left.

The words "Taping Socks" should be on the outer side of the foot.

The taping bands are divided into two function areas. The heel band goes on the outside, from the base of the Achilles tendon, under the outer ankle to the foot sole.

The upper edge of the taping band lightly touches the outer ankle ①.

The front taping band goes from the foot sole at a 45 degree angle to the ball of the foot at the back of the foot. The front edge lies behind the little toes ②.

On the inner side the band runs from the start of the Achilles tendon under the inner ankle towards the foot sole. Here the upper edge of the taping band lightly touches the inner ankle ③. Pulling the taping band stabilises the heel in the correct position and raises the inner longitudinal arch ④.

The front taping band on the inner side stabilises the first metatarsal bone towards the ground ⑤. Please use the two finger technique to check the position. Check that the lower leg is at a 90° angle to the foot by putting your index and middle fingers in the instep ⑥. If the back edge of the taping band is below your middle finger, it is correctly positioned ⑦.

Care instructions

Wash before use • wash inside out • do not use fabric softener

TAPE LOOP instructions for use

The rubber side of the taping band is applied to the skin.

Hallux valgus: Place the smaller loop around the big toe so that the bands cross immediately on the side behind the ball of the foot **1**. Make sure that the band does not cut in between the big toe and second toe.

Hammer toe: Place the smaller loop around the second toe or the second and third toes so that the bands cross beneath the ball of the foot **2**. Ensure that the loop of the band is positioned as far back on the toe(s) as possible to achieve the desired correction **3**.

Pull the big loop up behind the heel so that it cannot slide down. The taping band should be lightly tensioned and pull the big toe inwards / the small toe(s) out of their misaligned position and towards the ground.

Important: For optimal correction, you must ensure that the loops are sitting in the exact right position. If the correction tape cuts into the skin, use the next size up.

Care instructions

Use the COMPRESSION wash bag • do not use fabric softener

Use

The best results are achieved from COMPRESSION TAPE SOX and TAPE LOOPS if worn in shoes that are not too tight, with flat heels (< 4 cm) and enough room for the toes. The shoes must not be too short or too small, as this pushes the toes out of the correct position, which can cause bunions, as well as hammer/claw toes.

It is important that you slowly increase the length of time that you wear the corrective socks or tapes, or wear them first overnight, so that the toe joints can gradually get used to the new mobility. In the course of 1–2 weeks you can increase the length of time you wear COMPRESSION TAPE SOX or TAPE LOOPS until it is for the whole day.

Wearing time = treatment time

You can wear taping socks inside slippers, as well as outdoors and sports shoes. They can be worn for long times. As tabi socks they are suitable for wearing with sandals or without shoes.

The TAPE LOOPS can be worn with no limitations under normal socks, barefoot in summer shoes, or in sandals.

The familiarisation period

It's important that during the familiarisation period you also strengthen your muscle structure. There are various exercises that help to relax and strengthen your muscles.

An example for bunions:

Exercise to strengthen the big toe abductor

- Put your foot on the ground.
- Put more weight on the outer edge, particularly the outer ball than the inner side
- Brace the big toe inwards against the floor
- Hold the tension for 10 seconds
- Relax the foot
- Repeat the exercise 3–4 times

Try to do the movement in a controlled way, consciously using the muscles. Often the nerve impulse (conscious action) has been lost. So first it needs to be controlled, then strengthened as above.

Your physiotherapist will give you more exercises for your legs, feet and joints or you can find them on the (German-language) website www.my-medibook.de/uebungen-nach-diagnose.html (by Dr. Jens Wippert)

Simply scan
the QR code
for exercises.

Storage / Durability / Warranty

Store at room temperature and in a dry place. Avoid heat and moisture. Protect from sunlight.

The warranty period for COMPRESSANA TAPE SOX and LOOPS is 6 months starting from first use, as taping socks and tapes are subject to aging caused by material and wearing.

The warranty/guarantee also extends to the use of COMPRESSANA accessories.

In case of item damage please contact your specialist retailer or dispensing chemist. Do not attempt to mend the taping socks yourself as this could affect the therapeutical effect.

Risks and side effects

Observe the instructions of the medical personnel. Taping socks can cause skin necroses, oedems and stress damage on peripheral nerves, in particular with inappropriate handling. If pain or negative changes (skin provoking) during application should arise, refrain from wearing TAPE SOX and LOOPS and immediately contact your doctor.

In case of sensory and/or wound healing disturbances (e.g. diabetic diseases) implicating a skin lesion, badly healing impairments can develop unnoticed. Patients of this group require special control as well as a specific instruction about particular risks and precaution possibilities.

If there are incompatibilities materials and substance of content please consult your doctor, pharmacist or specialist retailer before use.

The manufacturer is not responsible for damage and injuries, which developed due to inappropriate handling, misuse or wrong handling due to non-consideration of our references or use of accessories not originating from COMPRESSANA.

Every foot is individual

Choosing the right size / correction strength and achieving a perfect fit is critical. The size / sizing information or correction strength is embossed on the tapes. If in doubt or in the event of any problems, please consult your specialist retailer and let them check the positioning and fit.

For hygienic reasons, this product is intended for one person only.

Disposal

No special disposal criteria exist. COMPRESSANA TAPE SOX and TAPE LOOPS can be disposed over the domestic refuse, the packaging material with recycled paper.

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COMPRESSANA GmbH