FOR COMPRESSION THERAPY



ASSORTMENT WITH SYSTEM

Top quality in terms of medical effect, optimum wearing comfort and fashionable aspects go hand in hand at COMPRESSANA – for healthy and beautiful legs!

In addition to compression care, COMPRESSANA also offers individual flat-knitted solutions and custom-made products as well as specially developed compression systems and bandage aids that are used for venous and lymphologic disorders such as venous leg ulcers, lymphoedema and lipoedema.

		NOVA	LEGANZA	DIAMONDS	SENSITIV	COTTON	LINEA COTTON SPORT COTTON	FREE	FORTE	ULCER SYSTEM/ INTRA	COTTON pro	FORTE pro	FORTE stabil
white			S/M	S	S/M								
racing gr	ey		C/M		C/M	S/M			C/M				
sandel		S/M	S/M S/M	S S	S/M S/M	S/M	S/M	S/M	S/M S/M	S/M	м	M	M
silk pearl		S/IVI S/M	S/IVI S/M	5 S	S/M	5/11	5/11	5/11	5/11	5/11	IVI	IVI	IVI
choco		3/11	S/M	S	S/M	S/M						М	M
graphi			S/M	S	S/M	S/M			S/M			M	M
indigo			S/M	S	S/M	S/M			S/M			M	M
black		S/M	S/M	S	S/M	S/M	S/M	S/M	S/M			M	M
Gradient Co	louro												
rising ind			S/M										
rising gra			S/M										
rionig gru	princ		0/11										
Trend Cold	ours		S/M	S	S/M								
					Der	und-kni [.]	Had					at Locate	bd
					ROL						E E	ат-конце	
	CCL 1	180	160	168				170	130	183 ³		at-knitte 135	
	CCL 1 CCL 2	180 280	160 260	168 268	195 295	156 256	150 250	170 270	130 230	183 ³	155 255	135 235	125 225
closed toe					195	156	150			183 ³	155	135	125
closed toe	CCL 2				195	156	150		230	183 ³	155 255	135 235	125 225
closed toe	CCL 2 CCL 3				195	156	150		230	183 ³	155 255	135 235 335	125 225 325
closed toe open	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2	280	260		195 295	156	150 250	270	230 330	183 ³	155 255 355	135 235 335 435	125 225 325 425
	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2 CCL 3	280 181	260 161		195 295 	156	150 250 	270 170	230 330 130	183 ³	155 255 355 155	135 235 335 435 135 235 335	125 225 325 425 125 225 325
open	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2	280 181	260 161		195 295 	156	150 250 	270 170 270	230 330 130 230		155 255 355 155 255 355	135 235 335 435 135 235 335 435	125 225 325 425 125 225
open	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2 CCL 3	280 181	260 161		195 295 	156	150 250 	270 170 270 A-d	230 330 130 230		155 255 355 155 255 355 Standard sizes in A-d	135 235 335 435 135 235 335 435	125 225 325 425 125 225 325
open toe Style	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2 CCL 3 CCL 3 CCL 4	280 181	260 161		195 295 	156 256	150 250 150 250	270 170 270 A-d A-g/H	230 330 130 230	3333	155 255 355 155 255 355	135 235 335 435 135 235 335 435	125 225 325 425 125 225 325
open toe Style Borders	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2 CCL 3 CCL 3 CCL 4	280 181	260 161		195 295 190 290	156 256	150 250 150 250 	270 170 270 A-d	230 330 130 230 330	333 ³ A-d	155 255 355 155 255 355 Standard sizes in A-d	135 235 335 435 135 235 335 435	125 225 325 425 125 225 325
open toe Style	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2 CCL 3 CCL 3 CCL 4	280 181	260 161 261	268	195 295 190 290 	156 256	150 250 150 250 Style S/M	270 170 270 A-d A-g/H 2 A-G -	230 330 130 230 330 5/M	3333	155 255 355 155 255 355 Standard sizes in A-d	135 235 335 435 135 235 335 435 Standard sizes in A-d CCL 1+2	125 225 325 425 125 225 325
open toe Style Borders Tricot bor	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2 CCL 3 CCL 4 CCL 3 CCL 4	280 181 281	260 161	268 	195 295 190 290	156 256 	150 250 150 250 	270 170 270 A-d A-g/H 2 A-G	230 330 130 230 330	333 ³ A-d	155 255 355 155 255 355 Standard sizes in A-d CCL 1+2	135 235 335 435 135 235 335 435	125 225 325 425 125 225 325 425
open toe Style Borders Tricot bor Lace bai	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2 CCL 3 CCL 4 CCL 3 CCL 4	280 181 281 - S/M	260 161 261 - S/M	268 	195 295 190 290 290 S/M S/M	156 256 	150 250 150 250 50 Style S/M S/M	270 170 270 A-d A-g/H A-G - S/M	230 330 130 230 330 5/M S/M	333 ³ A-d	155 255 355 155 255 355 51 Standard sizes in A-d CCL 1+2 - M	135 235 335 435 135 235 335 435 5tandard sizes in A-d CCL 1+2 - M	125 225 325 425 225 325 425 - M
open toe Style Borders Tricot bor Lace ba Softcontact	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2 CCL 3 CCL 4 CCL 3 CCL 4 S der nd band dots	280 181 281 - S/M -	260 161 261 - S/M S/M	268	195 295 190 290 290 5/M S/M	156 256 	150 250 150 250 5/M S/M S/M	270 170 270 A-d A-g/H A-g/H A-G - S/M S/M	230 330 130 230 330 5/M S/M S/M	333 ³ A-d	155 255 355 155 255 355 Standard sizes in A-d CCL 1+2	135 235 335 435 135 235 335 435 5tandard sizes in A-d CCL 1+2 CCL 1+2	125 225 325 425 125 225 325 425 425 - M M
open toe Style Border Tricot bor Lace bat Softcontact Silicone bar	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2 CCL 3 CCL 4 CCL 3 CCL 4 CCL 3 CCL 4 S der nd band dots / pattern	280 181 281 - S/M - S/M	260 161 261 - S/M S/M S/M	268	195 295 190 290 290 S/M S/M S/M S/M	156 256 	150 250 150 250 Style S/M S/M S/M S/M	270 170 270 270 A-d A-g/H A-G - S/M S/M	230 330 130 230 330 330 S/M S/M S/M S/M	333 ³ A-d	155 255 355 255 355 355 355 355 355 CL 1+2 - M M M	135 235 335 435 135 235 335 335 335 335 335 435 CCL 1+2 CCL 1+2 - M M M	125 225 325 425 125 225 325 425 425 - M M M
open toe Style Borders Tricot bor Lace ba Softcontact Silicone bar Silicone bar Silicone dots w	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2 CCL 3 CCL 4 CCL 3 CCL 4 der nd band dots / pattern s	280 181 281 - S/M - S/M	260 161 261 - S/M S/M S/M	268	195 295 190 290 290 S/M S/M S/M S/M	156 256 	150 250 150 250 Style S/M S/M S/M S/M	270 170 270 A-d A-g/H A-g/H S/M S/M S/M	230 330 130 230 330 330 S/M S/M S/M S/M	333 ³ A-d	155 255 355 255 355 355 355 355 355 CL 1+2 - M M M	135 235 335 435 135 235 335 335 335 335 335 435 CCL 1+2 CCL 1+2 - M M M	125 225 325 425 125 225 325 425 425 - M M M
open toe Style Borders Tricot bor Lace ban Softcontact Silicone bar Silicone dots w Borders	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2 CCL 3 CCL 3 CCL 4 CCL 3 CCL 4 der nd band dots / pattern s otos (3.5cm)	280 181 281 - S/M - S/M S/M	260 161 261 	268	195 295 190 290 5/M S/M S/M S/M S/M	156 256 	150 250 150 250 S/M S/M S/M S/M S/M S/M	270 170 270 A-d A-g/H A-g/H S/M S/M S/M S/M S/M A-D	230 330 130 230 330 S/M S/M S/M S/M	333 ³ A-d	155 255 355 255 355 355 Standard sizes m A-d ic A-d ic A-d ic A-d M M M M M M	135 235 335 435 235 335 435 Standard sizes in A-d icCL 1+2 - M M M M M M	125 225 325 425 225 325 425 - M M M M M
open toe Style Borders Tricot bor Lace ban Softcontact Silicone band d Silicone band d	CCL 2 CCL 3 CCL 4 CCL 1 CCL 2 CCL 3 CCL 3 CCL 3 CCL 4 der nd band dots / pattern S ots (3.5cm) nd (3.5cm)	280 181 281 - S/M - S/M S/M	260 161 261 	268 	195 295 190 290 290 5/M 5/M 5/M 5/M 5/M	156 256 	150 250 150 250 Style S/M S/M S/M S/M S/M S/M	270 170 270 270 A-d A-g/H 3 A-G - S/M S/M S/M S/M S/M S/M	230 330 130 230 330 5/M S/M S/M S/M S/M S/M		155 255 355 155 255 355 Standard sizes in A-d CCL 1+2 M M M M M S ¹ /M	135 235 335 435 235 335 435 Standard sizes in A-d CCL 1+2 - M M M M M S ¹ /M	125 225 325 425 225 325 425 425 - M M M M M

AN OVERVIEW OF THE RANGE

PRODUCTS FOR COMPRESSION THERAPY - round knitted	4
COMPRESSANA LEGANZA – incomparably attractive	4
COMPRESSANA DIAMONDS – with diamond pattern	7
COMPRESSANA NOVA – uncompromisingly transparent	8
COMPRESSANA SENSITIV - incredibly smooth	10
COMPRESSANA COTTON – consequently natural	12
COMPRESSANA LINEA COTTON - rib look	14
COMPRESSANA SPORT COTTON – sporty and modern	15
COMPRESSANA FREE – easy and comfortable	16
COMPRESSANA FORTE – reassuringly effective	18
COMPRESSANA ULCER SYSTEM	20
PRODUCTS FOR COMPRESSION THERAPY – flat knitted	22
COMPRESSANA FORTE pro	24
COMPRESSANA FORTE stabil	25
COMPRESSANA COTTON pro	26
COMPRESSANA SEGMENTS	27
	21
COMPRESSANA SLEEVES, HAND + TOE PARTS	28

COMPRESSANA COMPRESSION PADS	29
INTRA BASE LAYER STOCKINGS	30
MEDICAL ADAPTIVE COMPRESSION SYSTEMS (MAC)	32
COMPRESSION WRAP ACS Light	34 35
BANDAGEHILFEN SOFTCOMPRESS	36
SCAR COMPRESSION PATCH ScarExPress	37
MASTER Donning and doffing aid SNAPPER Donning Aid EASY Donning Aid SPECIAL GLOVES FIXATION	38 39 40 41 42
LYMPHATIC CUSHION	43
SKIN CARE	44
ACCESSORIES FOR CARE, WASHING AND STORAGE INSHOE-SOCKS AND STEPS	47 49

COMPRESSANA LEGANZA

incomparably attractive



TENDER TO LEGS ATTRACTIVE IN APPEARANCE

- exciting, highly sensitive
- extremely fine and silky shining
- even, smooth surface structure

By the use of high tech yarns and due to an optimized manufacturing and finishing process, an unequalled homogeneous and fine mesh surface is obtained. COMPRESSANA Leganza convinces with a favourable smooth surface!

The very fine texture causes a pleasant feeling and a marvellous softness on the skin and makes COMPRESSANA Leganza smooth and soft.

COMPRESSANA Leganza offers a superior appearance due to its innovative material properties. Compression yarns in its natural state facilitate a discreet transparency and accentuate thereby the natural skin tone.

The shape, which is specially adapted to the anatomy, and the excellent stretch properties ensure optimum wearing comfort in every situation. Furthermore COMPRESSANA Leganza convinces with a favourable accentuation of the legs and reveals a feminine silhouette.


Style

A-d

CCL 1 light compression 2.4 - 2.8 kPa (18 - 21 mmHg) Art.No. 160 (closed toe) Art.No. 161 (open toe)

CCL 2 moderate compression 3.1 - 4.3 kPa (23 - 32 mmHg) Art.No. 260 (closed toe) Art.No. 261 (open toe)







NEW: Softcontact band with microdots for best grip with low pressure load

Colours			
white	sandel	silk	pearl
choco	graphit	indigo	black

COLOURS Compression valid from 02/2022							
lilac	pacific		grape				
Gradient Colours valid from 02/2022							

rising indigo

rising graphit

Extravagant colour gradient decreasing towards the heart for thigh-high stockings and tights

COMPRESSANA LEGANZA & DIAMONDS Perfect fit - Superior wearing properties




Panty

- · anatomically shaped panty
- different elasticity and stretch zones
- thin meshes and non-bulky fabric
- flat waistband, convenient and comfortable



Foot

- completely reinforced sole without restriction of the longitudinal elasticity
- spacious, non-constricting toe with soft stretch
- soft and particular fine seams at the toe



Knitting

- natural yarns for feminine gloss and delicate appearance
- softness, transparency and unusual homogeneous and fine mesh surface by innovative yarn and knitting technology
- knitted on most modern knitting machines



Thigh zone

- flexible, high elastic thigh zone with stretch reserves
- feminine silicone top bands in Italian design, sewn without overlapping



Knee and calf zone

- soft, elastic, non-constricting in the hollow of the knee
- · below-knee stockings with seamless soft-cuff

Convincing in application and handling

- streak-free and without colour transitions
- gently coloured and elaborately refined
- climate, temperature and moisture regulating characteristics
- easy-care, suitable for machine wash and tumble dryer


DIAMONDS

with diamond pattern

DIAMONDS ARE A GIRL'S BEST FRIEND

COMPRESSANA DIAMONDS with knitted diamond pattern and fine silky mesh shows your egs confidently. The modern way of compression therapy: Invisible effect – visibly attractive! COMPRESSANA DIAMONDS is available in standard sizes.









Styl	е					
P	A-d A-d short A-d extra short	Below-k	nee stocki	ings		
P	A-d A-d short A-d extra short	 Silicone 	nee stocki band (3.50 ntact band (cm)		
P	A-g/H A-g/H short A-g/H extra short		gh stockin Itact band	igs		
Ł	A-t A-t short A-t extra short	Tights asymmetric panties, Body extends standard size approx. 80-120 cm				
Ł	A-t/G A-t/G short A-t/G extra short	Maternit	y tights			
COLOURS Compression valid from 02/2022						
lil	ac pacific	stone	grape	sunrise		



NOVA

uncompromisingly transparent



TRANSPARENCY IN COMPLETION

- particularly transparent
- airy and fine knitting
- ideal in summer

durable

·BASIC

attractiv

Inconspicuous appearance, air permebility, transparency and perfect fit makes COMPRESSANA NOVA unique among compression stockings. Our customers also think so:

I would like to say a big praise for NOVA at this point. So far - and I have tried quite a few models over the years - these are the first stockings that are really comfortable to wear, that last all day without me having to pull them up. I can wear them even in warm temperatures without feeling uncomfortable. They are also much less susceptible to running. Thank you so much for bringing this model to the market, you have done a great service to many people and made compression therapy HIGHLY more comfortable!

As an addendum, I would like to say that NOVA stockings are also the best in terms of fashion. I'm 39 and I certainly don't feel like walking around in shiny or totally opaque granny stockings. I wear NOVA normally with skirts and it is not at all noticeable that they are compression stockings.

Many greetings M.B.



attractive lace-band




anatomical perfect toe with flat seam



CCL 2 moderate compression 3,1 - 4,3 kPa (23 - 32 mmHg) Art.No. 280 (closed toe) Art.No. 281 (open toe)





Colours		
silk	pearl	black

SENSITIV

incredibly smooth

-



VELVET AND UNBELIEVABLY GENTLE

- incomparably soft and gentle
- especially comfortable to wear
- temperature and moisture regulating, breathable

COMPRESSANA SENSITIV coats the leg with a finely woven and soft knitted fabric. COMPRESSANA SENSITIV is easy to put on and grants maximum wearing comfort. The compression spreads gently and evenly. Annoying pinching and feeling of pressure does not exist any longer. Instead, COMPRESSANA SENSITIV provides a supple feeling on the skin thanks to its soft and smooth properties – you immediately forget that you are wearing a compression stocking. The high-quality, double covered yarns flatter the skin and offer best skin compatibility, skin protection and long durability.


The stockings are breathable and moisture regulating. The even and concealing fabric offers both ladies and men a visually attractive appearance.



CCL Art.No. 295 (closed toe) Art.No. 290 (open toe)



Styl	e	
P	A-d A-d short A-d extra short	Below-knee stockings
P	A-d/H A-d/H short A-d/H extra short	Below-knee stockings • Silicone band (3.5 cm) • Softcontact band (3.5 cm)
P	A-g A-g short A-g extra short	Thigh-high stockings
		Thigh-high stockings • Silicone band (dots) • Lace band • Silicone dots w/ pattern • Softcontact band
P	A-g/T A-g/T short	Thigh-high stockings with hip attachment
Þ	A-t A-t short A-t extra short	Tights Highly elastic pantyhose with adjustable waistband Body extend approx. 80-140 cm
P	B-t B-t short	Leggins in CCL 1 (for combined use)



COLOURS Compression valid from 02/2022						
lilac	pacific	stone	grape	sunrise		

COMPRESSANA

consequently natural




THE NATURAL SOLUTION FOR LEG PROBLEMS

- soft, smooth cotton inside
- air-permeable and very skin friendly
- hardwearing by abrasionresistant exterior

COMPRESSANA COTTON is the skin-friendly option for compression stockings. The high cotton content ensures a particularly pleasant and natural wearing comfort. COMPRESSANA COTTON is suitable for people who suffer from dry and itching skin or who prefer natural materials.

The special feature of COMPRESSANA COTTON is its processing. Due to a special technology, the cotton is only on the inside of the stocking. This results in an almost 100% coverage of the legs with natural fibres, while the wear-resistant synthetic fibres protect the outside of the stockings against fast abrasion.

COMPRESSANA COTTON is the ideal compression stocking for those who want a natural feeling on the skin.





open toe

black



COMPRESSANA LINEA COTTON

COMPRESSANA LINEA COTTON with ribbed structure provides below-knee stockings with a modern and sporty appearance. Inconspicuously, so that nobody can recognize by the mesh that you are wearing compression stockings.


COMPRESSANA LINEA COTTON is also available with grip top for double secure grip. The silicone beaded band is slip-resistant and suitable for highest demands. COMPRESSANA LINEA COTTON is ideal for moving, active legs due to the unique and skin-friendly characteristics and its high cotton content.

COMPRESSANA LINEA COTTON is advantageous to other compression stockings. The graduated compression and the efficient depth effect increases performance and persistence.





closed toe



Colours				
silk	choco	graphit	indigo	black

SPORT COTTON

Performance and endurance

The wearing comfort and the outstanding properties of LINEA COTTON in racing style for all ambitious athletes and sportsmen.

- · better blood circulation and oxygen supply
- activation and stabilisation of the musculature
- faster recovery
- reducing of the risk of injury



moderate compression 3.1 - 4.3 kPa (23 - 32 mmHg) Art.No. 256



 A-d
 Below-knee stockings

 A-d short
 with ribbed structure

 A-d/H
 with ribbed structure

 A-d/H
 Below-knee stockings

 A-d/H
 silicone band (3.5 cm)

 Softcontact band (3.5 cm)
 Softcontact band (3.5 cm)

Colours





COMPRESSANA FREE light & comfortable

FOR MORE FREEDOM IN THERAPY


- unsurpassed ease of donning and doffing
- more mobility
- pleasant and comfortable
- better hold, less slipping

Innovative knitting technology with integrated functional zones, sliding yarn technology and controlled elasticity for a new feeling of compression. Due to its diverse properties and benefits, COMPRESSANA FREE is suitable for a very wide range of indications for men and women of any age and constitution and offers noticeable therapy relief even for beginners or people with physical disabilities/restrictions.

In addition to its look suitable for everyday use, COMPRESSANA FREE with its soft, skin-friendly and breathable yarns ensures optimal wearing comfort, more freedom of movement and uncomplicated handling, making it easy to put on and take off and promoting mobility.

COMPRESSANA FREE enables noticeably easier compression therapy with maximum comfort – when putting on and wearing.





nnovative functional zones

- \cdot for an easy & pleasant everyday life
- for a better wearing feeling & noticeable therapy efficiency




FORTE reassuringly effective

STRONG AND STURDY IN EVERY CONDITION OF LIFE

- durable and long lasting
- easy donning
- soft and smooth
- outstanding elasticity
- extraordinary in the knee and thigh-part

The stocking for everybody for whom function counts most. With COMPRESSANA FORTE you uncompromisingly decide for a stocking with exceptional effectiveness with which you walk through thick and thin. Strong, covered varns create firm vet comfortable compression stockings. With strong working pressure, COMPRESSANA FORTE is also suitable for men. The wide-meshed and soft knit promotes air circulation, transports moisture and stimulates blood circulation.

Another advantage: COMPRESSANA FORTE is particularly easy to put on, as the robust, wear-resistant material is easy to handle.





open toe

strong compression **4**.5 - 6.1 kPa (34 - 46 mmHq)

black

open toe

COMPRESSANA ULCER SYSTEM

simply practical

WITH COMPRESSION TO SUCCESS

Medical 2-component compression stocking system for the treatment of venous leg ulcer In the case of venous leg ulcers, after the decongestion phase of the oedema, the treatment should be switched to two-layer ulcer compression systems for long-term therapy, if possible. The COMPRESSANA ULCER SYSTEM consists of an overstocking and two understockings that can be changed. When worn in combination, the compression pressure of ULCER FORTE and ULCER INTRA add up to a resting pressure of compression class 3. The stocking components are optimally matched to each other and easy to use: The understocking with light compression is used to fix wound dressings and pads. It can be worn on the affected leg day and night for permanent compression. The overstocking with medium compression is applied over the understocking during the day to increase compression.

Compression in combination with movement ensures treatment success: COMPRESSANA ULCER SYSTEM offers non-slip wearing comfort, full freedom of movement and mobility. Normal sports clothing and shoes can be worn. Even after the venous leg ulcer has healed, medical compression stockings with high working pressure and depth effect remain indispensable. The COMPRESSANA ULCER SYSTEM can also be used for recurrence prophylaxis.

strong Compression 4.5 - 6.1 kPa (34 - 46 mmHg) Art.No. 333

ULCER SYSTEM

1 Overstocking ULCER FORTE CCL 2 with open toe + 2 Underliner ULCER INTRA CCL 1 with closed toe

Colours	
silk	



ULCER INTRA (Additional delivery set) 2 Underliner ULCER INTRA CCL 1 with closed

toe

(Made-to-measure: 1 Underliner)

Style



Below-knee stockings



Overstocking with open toe and Underliner



Underliner: Heel / instep marking for exact positioning



HIGHER WEARING COMFORT AND PERFECT FIT STITCH BY STITCH

THERAPEUTIC FLAT KNITTED COMPRESSION GARMENTS

Compression therapy focuses on venous diseases (phlebology). For this treatment, round knitted compression stockings are used in the predominant compression classes 1 and 2 (CCL 1 and CCL 2). Here primarily ready to wear sizes for anatomically normal legs are applied.

Flat knitted compression products up to CCL 4 are almost only possible as custommade products.

Field of application:

- in lymphology for the medical ctreatment are of lymphoedema and lipoedema, as well as in mixed forms e.g. lymphoedema combined with venous disease
- scar treatment, especially during therapy of deep burns and scalds
- treatment of anatomically complex and problematic body proportions that cannot be treated with circular-knitted stockings
- prosthetic garments for shaping oedematous soft tissue as well as for the improvement of wound healing and pain reduction after leg amputations

The advantages of flat knitted compression stockings:

- optimum, precise fit
- stronger effect through high working pressure during movement
- better grip and less pressure points and rubbing
- high stability ensures fit without folds and constrictions
- enhances subcutaneous microcirculation through massage effect



Stitch by stitch


Flat knitted compression stockings are manufactured stitch by stitch in a complex process. Only the flat knitting production process allows the number of stitches to be altered within the stitch course. Also, modifications can be made including but not limited to instep stitching, anatomical flexes, homogenous endings and slanted tops/toes. After completion of the knitted fabric, both endings are sewn together or better the stitches on both edges of the seam are linked with a single thread. This ensures optimum flexibility of the joints under compression.

Flat knitted compression products have less tension and are dimensionally stable. They do not lay into skin folds and prevent undesirable constrictions. The manufactured compression stockings or compression model the anatomy exactly even in case of vast proportional differences.



FORTE pro

COMPRESSANA FORTE pro has a particularly supple knitted fabric that is characterised by its high wearing comfort and uniform appearance. It allows full freedom of movement and is therefore ideal for the permanent/long-term care of unproblematic or fully treated oedema. Due to the relatively flat, fine-meshed knitted fabric, FORTE pro is considered to be rather comfortable, light and pleasant compared to stable variants. Thanks to the high UV absorption, the regeneration of young skin is supported. This property is particularly advantageous in scar and burn therapy.




FORTE stabil

Stable and strong: these are the properties of COMPRESSANA FORTE stablil. The high material thickness and strength ensure a stable wall effect and an unsurpassed depth effect. The air-permeable, porous, three-dimensional knitted surface develops a strong massage effect during movement.

The material thickness and stiffness of the knitted fabric ensures an even distribution of the contact pressure in case of calibre jumps e.g. before/after joints, skin furrows or skin flaps.








COMPRESSANA COTTON pro

COMPRESSANA COTTON pro wears naturally and is very skin-friendly thanks to its high cotton content. Increased moisture absorption and excellent temperature regulation make it a reliable everyday companion. It allows full freedom of movement with effective working pressure and depth. The soft and slippery material makes COMPRESSANA COTTON pro particularly recommended for **sensitive skin**, skin irritations or for patients who prefer cotton on the skin.

light compression 2.4 - 2.8 kPa (18 - 21 mmHg) Art.No. 155

moderate compression 3.1 - 4.3 kPa (23 - 32 mmHg) Art.No. 255


strong compression 4.5 - 6.1 kPa (34 - 46 mmHg) Art.No. 355



flat-knitted with natural cotton

• skin-friendly


 effective in penetration and working pressure



SEGMENTS AND APPLICATION AREAS

of flat knitted products

Combinations possible within the framework of two-piece or multi-piece compression fittings for arm/hand and leg/foot





COMPRESSANA SLEEVES, HAND- & TOE PARTS







ccL 3 strong compression 4.5 - 6.1 kPa (34 - 45 mmHg)



COMPRESSANA COMPRESSION PADS

In compression therapy, pressure pads are used on the palms/back of the hands or on legs in the ankle pits to pad concave areas, to create the most even pressure possible or to increase pressure in a targeted manner. Pads are mainly used in oedema therapy to promote microcirculation, increase lymphatic drainage, resorb oedema, remove deposits and loosen fibrosis. Also in scar treatment, compression promotes tissue healing, reduction of

scar thickness and reduction of itching and redness of the skin. In addition to silicone sheets and plates, pads are also used to increase or optimise pressure.

COMPRESSANA pads offer individual fitting options with the best possible fitting, optimal wearing comfort and skin tolerance!



Possible applications under compression stockings, segments, bandages and dressings:

- Loosely inserted (individual placement by the customer) or fixation with COMPRESSANA ULCER SYSTEM or COMPRESSANA INTRA.
- Sewn in (glued or sewn): Is glued into the stocking/bandage with silicone glue at the chosen location after being marked by the customer and then covered with lining fabric so that the location remains elastic.
- Slipped into pockets made of lining fabric: Pocket is sewn on the inside of the sock/segment after marking by the customer and pad can be slipped in.



INTRA toesock with forefoot support Art. 0090



Half toe sock, open toe Cotton/Silver Art. 0095



Base layer sock with Silver Art. 0094



Base layer below-knee stocking with Silver Art. 0096/0086



Base layer below-knee stocking with Cotton Art. 6082



Base layer below-knee stocking COMPRESSANA ULCER INTRA Art. 0183



Base layer thigh-high stocking with Cotton Art. 6082

COMPRESSANA INTRA

Base layer and fixing system for modular wound and skin treatment

Due to its soft knitted fabric, INTRA is particularly comfortable to wear and is suitable as an understocking under medical compression stockings and systems. The stockings have no compression and are therefore particularly easy and gentle to put on and take off. INTRA offers permanent skin protection and has an additional antibacterial effect thanks to the knitted-in silver thread. It is ideal for fixing special pads and for protecting wound dressings when putting them on and taking them off. The liner protects sensitive skin and, if necessary, fixes the wound dressing during donning and doffing. The understocking is perfect for use in combination with Medical Adaptive Compression Systems.

Scope of application:

- Slipping-free fixation of truss pads, inlays to protect against pressure and pads for the treatment of medically and anatomically complex situations on the foot, ankle or lower leg.
- Toe parts (toe support or toe protection) with cotton to protect toes and interdigital spaces.
- Due to high flexibility suitable for wearing in preference to socks/below-knee stockings in presence of tenderness, deformation, adhesion, large leg, foot circumferences and excessive foot sweat, foot odour, skin or fungal disease. Can also be used as an understocking.

Product	Colour	Size	Art.No.	RRP*
INTRA Toe sock, closed toe with cotton with forefoot support (e.g. in case of hallux valgus)	white	One-Size	0090	12,90 €
INTRA Half toe sock, open toe	white	One-Size	0095	7,50€
INTRA Base layer sock, SkinProtect®Silver X-Static	silver	One-Size	0094	18,50 €
INTRA Base layer stocking, SkinProtect®Silver X-Static	silver	One-Size	0096	22,-€
INTRA Base layer stocking extra wide, SkinProtect®Silver X-Static	silver	Calf circumference > 50cm and large foot sizes	0086	23,-€
INTRA base layer below-knee compression stocking CCL1	silk	I – VI	0183	49,-€
Liner base layer below-knee stocking with cotton	black	1 (calf circumference 29-48 cm) 2 (calf circumference 48-69 cm)	6082	29,90 €
Liner base layer Thigh-high stocking with cotton	black	1 (tigh circumference 50-75 cm) 2 (tigh circumference 76-99 cm)	6082	39,90 €
* approx recommanded retail price				

* approx. recommended retail price

MEDICAL ADAPTIVE COMPRESSION SYSTEMS (MAC)

- easy and safe to use
- less time and effort required
- faster treatment success
- more freedom of movement and independence
- more self-determination and personal responsibility
- better quality of life

COMPRESSION PRODUCTS FOR DECONGESTION


Self-bandaging made easy

People suffering from oedema and its effects need decongestive therapy to reduce swelling and relieve discomfort. Above all, compression bandaging ideally ensures the success achieved by therapists. In contrast to conventional wrapping bandages, Medical Adaptive Compression systems (MAC) and bandage aids can be applied more effectively and also independently by the patient, relatives or carers.

Innovative aids for self-bandaging are offered for a noticeable increase in quality of life. Compression systems are practical to use especially during the decongestion phase or if needed further during the maintenance phase. They can be worn day and night. They enable the patient to make individual adjustments while wearing them. Furthermore, depending on the situation, they can be ideally adapted to fluctuating body sizes. The success of the decongestion treatment can be further promoted by using the bandages for padding. Sensitive parts of the body are also protected from excessive shear and friction forces.

Medical Adaptive Compression Systems can also complement the entire therapy from decongestion to maintenance on a patient-specific basis in cases of leg ulcers, secondary diseases, intolerances and in cases that make compression therapy difficult.

The complex manual decongestive Therapy (CDT) – a 2-phase treatment



Phase 1: Decongestion phase

The decongestion phase is aimed at the mobilisation of the oedema. In the first phase of this treatment, all the components of complex manual decongestive therapy should be used 1 – 2 times daily. The education, training and involvement of the patient in his / her treatment is of particular importance in this phase because the measures are time-consuming an require active patient participation.

Phase 2: Maintenance phase

In the maintenance phase, the components of complex manual decongestive therapy are applied according to diagnosis and complaints, to conserve and optimize the treatment success of phase 1. This phase may have to repeated, depending on the progression of the illness – in particular for concommitant conditions. The patient should also be actively involved in the second phase of treatment.



COMPRESSION WRAP

Thanks to its exceptionally short-stretch properties, the Compression Wrap has a very high working pressure and a low resting pressure. The double-sided Velcro fasteners and the helpful grip zones make it easy to put on and can be ideally adapted to fluctuating body circumferences. The special not-ches on the flaps enable a gapless closure.

- · Flexibly adjustable from 30 to 60 mmHg compression pressure
- \cdot Need for very high working pressure (e.g. venous leg ulcers, pronounced fibrosis)
- · Very pronounced lymphoedema (e.g. stage II III lymphoedema)
- \cdot Grip zones at the upper and lower ends facilitate application
- \cdot Colour can be individually changed from beige to black
- · Cost-effective and time-saving in contrast to wraparound bandages

The individual segments can be easily combined as needed.







ACS LIGHT

With short-stretch properties, ACS Light can be used efficiently and even more comfortably for decongestion. Easy to use due to incorporated fixation aid. Adhesive straps at the upper end ensure secure hold and mobility.

- · Individually adjustable up to 40 mmHg compression pressure
- · Also suitable for secondary conditions that require little pressure or when high compression pressure is not tolerated (e.g. pain sensitivity, sensibility, stimulus perception)
- · Conditions requiring more flexibility and elasticity
- · Temporary increase in pressure in addition to compression stocking in the maintenance phase
- \cdot Anatomical fit for maximum comfort
- · Particularly durable
- · Great time saving in contrast to wraparound bandages

The individual segments can be easily combined as needed.











SOFTCOMPRESS BANDAGE AIDS

Underpadding for MAK, practical in decongestion therapy and compression bandaging

Mit SoftCompress kann der Erfolg der Entstauungstherapie mit medizinischen adaptiven Kompressionssystemen deutlich gesteigert werden. Die Unterpolsterung kann zudem Druckstellen und Hautschäden vorbeugen und die Lebensqualität fördern.

Die Produkte zeichnen sich im Gegensatz zu häufig verwendeten Synthetikmaterialien durch eine besonders hautfreundliche Baumwoll-oberfläche aus. Diese ermöglicht einen angenehmen Tragekomfort für den Patienten. Die vertikalen Schaumstoffkanäle bilden eine einzigartige Oberflächenstruktur, die die Ödemreduzierung unterstützt.

Die komfortablen Bandagehilfen können unter medizinischen adaptiven Kompressionssystemen (ACS Light oder Compression Wrap), unter Kompressionsbestrumpfung oder unter Kompressionsverbänden getragen werden und sorgen für zusätzliche Kompression und Polsterung sowie eine ideale Druckverteilung.

SoftCompress complete sets

SoftCompress LymphAid and BandageAid offer a complete solution for lymphological compression bandaging. LymphAid comes with a ready-to-use bandage and is therefore also suitable for self-bandaging. The more economical BandageAid contains SoftCompress Compression Rolls.



SCAR COMPRESSION PATCH

COMPRESSION	
~	

Prevention, correction and reduction of keloids and hypertrophic scars in the size range of the compression element.


Partial compression through compression therapy is one of the recommended and most effective treatment methods in scar therapy.

ScarExPress contains an innovative compression element at its core, embedded in a skin-friendly and waterproof tape with strong adhesive power. Continuous pressure prevents disorderly and excessive connective tissue remodelling and the formation of proliferating and raised scars (hypertrophic scars, keloids). Pressure softens scar tissue, reduces scar thickness and redness, even in older scars.

www.narbenpflege.online

Product	Compression area	Size	Content	Art.No.	PZN	RRP.
ScarExPress Scar Compression Patch S	1,5 x 1 cm	S	4 pcs.	5401	14817383	49,-€
ScarExPress Scar Compression Patch M	3,2 x 1 cm	Μ	3 pcs.	5402	14817408	49,-€
ScarExPress Scar Compression Patch L	4,9 x 1 cm	L	2 pcs.	5403	14817414	49,-€

WITH PRESSURE AGAINST SCARS





ScarExPress Video

MASTER

Donning and doffing aid



Donning

In a sitting position, without bending down and without exerting force by pre-stretching the stocking material.



Doffing

Removing the stocking is facilitated by a retaining element that secures the stocking for removal and a handle that "extends" the arms. COMPRESSANA Master facilitates donning and doffing of stockings • Hassle free and easy to use while sitting down • Without strenuous bending • Without exertion of force • Makes starting compression therapy easier • Ensures that the stockings are donned and positioned correctly

Prior stretching of the stocking reduces friction and protects the skin. In addition, this protects the fabric of the stockings against damage. COMPRESSANA Master is a simple device suitable for all people in special life circumstances e.g. restricted mobility etc.

The COMPRESSANA Master makes starting compression therapy easier, additionally it ensures the stockings are donned and positioned correctly. Donning and doffing aids are recommended or prescribed by the doctor if the donning of the stockings without therapeutic appliances is too difficult.

For example in cases of illnesses or disabilities with considerable dysfunctions of the finger, hand and arm joints as well as for people with limited finger grip.



Product	Art.No.	RRP
MASTER	0560	89,-€
MASTER XL Ankle circumference cB > 26 cm	0565	89,-€
Surface disinfection wipe, without alcohol, pack of 10	0614	6,90 €

SN **PPER** Donning Aid

Putting on and taking off according to the Snapper concept is quick and requires little effort. This is ensured by the special construction and the gliding ability of the Snapper donning facilitator.

It is also advantageous to use stocking gloves with grippy palms, as recommended in the stocking manufacturers' instructions for use. These are available as accessories.

COMPRESSANA Snapper protects the skin from injuries and stocking stitches from damage

due to mechanical effects, e.g. friction, finger/ toe nails, rings or jewellery, when putting on and taking off, especially in combination with Top-Grip stocking gloves. The stockings thus retain their therapeutic effect for longer.

COMPRESSANA Snapper makes it easier to start compression therapy and ensures that stockings are put on correctly. COMPRESSANA Snapper is suitable for all people who prefer it lighter or who have special circumstances or life situations, e.g. limited mobility.



Video instructions at www.compressana.de or scan QR code with smartphone/tablet.



SNAPPER/MASTER live

www.youtube.de/ compressana

Product	Art.No.	RRP
SNAPPER donning aid for open and closed toe	0089	36,-€
Surface disinfection wipe, without alcohol, pack of 10	0614	6,90€





COMPRESSANA EASY – Donning Aid Leg & Donning Aid Arm

COMPRESSANA EASY, the simplest and cheapest donning aid and an alternative to expensive racks The particularly slippy material goes around the foot and leg, enables donning of stockings even with high pressure without any problems Equally suitable for open and closed toe. COMPRESSANA offers also a simplification for the donning of compression sleeves The particularly slippy material enables donning of sleeves even with high pressure without any problems. · easy to handle

sturdy

· perfect for mobile nursing services

Product	Art.No.	RRP
EASY Leg HMV-Nr. 02.40.01.3041	0098	12,-€
EASY Arm One size HMV-Nr. 02.40.01.3040	0120	19,50 €



COMPRESSANA Sliding Aid

The sliding sock is pulled over the forefoot; its slippery material makes it easy to pull the compression stocking onto the foot. The sliding sock is then simply removed through the toe opening. When the toe is closed, the sliding sock can be removed with the toe flap pointing back towards the heel.

Product	Art.No.	RRP
EASY Sliding Aid	0100	3,-€

We make it easy for you to get dressed and undressed: Video instructions at www.compressana.de or scan QR code with smartphone/tablet.

SPECIAL GLOVES

TIP The use of gloves makes it much easier to put on stockings and trousers and protects the stocking stitches from damage by fingernails or jewellery.





COMPRESSANA SUPER-GRIP

Extra fine, skin-friendly and breathable textile gloves Extra grip, perfect fit, seamless, better fingertip feeling.

Product	Art.No.	RRP
SUPER-GRIP Gr.II (S) until 8cm handwidth*	0112	6,-€
SUPER-GRIP Gr.III (M) until 9cm handwidth*	0112	6,-€
SUPER-GRIP Gr.IV (L) until 10cm handwidth*	0112	6,-€
SUPER-GRIP Gr.V (XL) until 11 cm handwidth*	0112	6,-€

* The size is determined by measuring the hand width at the finger baseline



COMPRESSANA TOP-GRIP

Breathable, sensitive and durable textile gloves made of nylon with nitrile coating – seamless, hardwearing.

Product	Art.No.	RRP
TOP-GRIP Gr.II (S) until 8 cm handwidth*	0108	7,-€
TOP-GRIP Gr.III (M) until 9 cm handwidth*	0108	7,-€
TOP-GRIP Gr.IV (L) until 10cm handwidth*	0108	7,-€
TOP-GRIP Gr.V (XL) until 11 cm handwidth*	0108	7,-€



COMPRESSANA GRIP

Special gloves made of natural rubber (latex) with non-slip palm, very good dexterity and cotton lined.

Product	Art.No.	RRP
GRIP Gr.II (S) until 7 cm handwidth*	0102	3,50 €
GRIP Gr.III (M) 7 to 8 cm handwidth*	0102	3,50 €
GRIP Gr.IV (L) 8 to 9 cm handwidth*	0102	3,50 €

FIXATION

The right fixation for every stocking

Integrated fixing aids, silicone-coated bands, internal bands, hip fixation, bodice, shoulder strap or bolero are available for fixing in the leg and arm area.



Silicone beaded band



Lace band



Silicone dots w/ pattern



Softcontact band

The attractive and the neutral silicone bands of the highest quality are breath-able, comfortable to wear and adapt perfectly to the leg. Choose according to your needs and preferences, by design and type of coating.



Hip and body attachment

Skin adhesive for effortless hold

The skin-friendly adhesive lotion, which can be removed without leaving any residue, gives stockings, bandages and other textiles a secure and long-lasting hold. Simply spread the adhesive

lotion on the skin, press the stocking or bandage against it - and you're done. After wearing, it can be removed from the skin within seconds with water or a damp cloth.



Produkt	Art.No.	RRP
Adhesive lotion Roll-On 60 ml	0126	13,-€
COMPRESSANA FIX Adhesive lotion Roll-On 60 ml, alcohol free	0116	13,50 €

LYMPHATIC CUSHION

Positioning and supporting aid for arm and hand

Comfort, hygiene and maximum effectivity - also during sleep and at work

The healing of lymphatic oedema and serious injuries of the arm and hand is supported by an elevated resting position. Swelling is reduced or even prevented Improves back flow of the blood to the heart as well as the removal of the lymph fluid. The COMPRESSANA lymphatic cushion is a valuable treatment aid for "complex physical decongestive therapy" (CPDT). It can also be used for pain treatment and wound management. Patients with lymphatic oedema, lip oedema or burns/scalding quite often suffer from bacterial infections (erysipel, abscesses, pus-producing infections and staphylococcus). These infections cause additional swelling, itching, redness of the skin and pain Therefore, the armrest surface of the COMPRESSANA lymphatic cushion is coated with a silver lining. The silver lining has a bacteriostatic effect, keeps the natural balance of the skin flora and can help to reduce irritating side effects of infections. Additionally, it protects the affected skin parts against contracting further infections.



Product	Art.No.	RRP
LYMPHATIC CUSHION Gr.I \leq 30 cm Size is determined by the circumference of the forearm	0550	179,95 €
LYMPHATIC CUSHION Gr.II > 30 cm Size is determined by the circumference of the forearm	0550	179,95 €





SKIN CARE

For care of skin under treatment with compression garments

Optimal care products for compression therapy

To ensure successful compression therapy for vascular diseases, oedemas etc., a healthy, elastic and resistant skin is of great importance. The development of the products was carried out in close cooperation with renowned dermatological cosmetics manufacturers. We exclusively use formulas with ingredients that have successfully proven themselves in practise for many years in accompanying therapy at clinics and practitioners. With the product line COMPRESSANA SKIN CARE COMPRESSANA offers care products especially adapted to the needs of compression therapy for prophylaxis and skin protection.

> For more information on Skin Care, please see our customer information (Art.No.: 0737).



Further information on Skin Care & Compression can be found in the download area atwww.compressana.de

COMPRESSANA ANTI-ITCHING SPRAY

The active ingredient polidocanol soothes the skin and can prevent itching for several hours. The spray provides moisture and has an antimicrobial effect against excessive bacteria and germ formation. Inflammation caused by scratching is prevented.

It can be used several times a day on affected skin areas or directly on support and compression stockings, is absorbed without residue and does not leave any stains. The spray can be used independently or in addition to the cream or lotion.

Ingredients: Aqua, Butylene Glycol, Laureth-9, Chlorhexidine Digluconate

Product	Art.No.	RRP
ANTI-ITCHING SPRAY Spray 150 ml	0586	15,-€
ANTI-ITCHING SPRAY Spray 50 ml	0580	10,-€

made of

plastic

RECYCLED

COMPRESSANA TONIC

- Vitalizes + strengthens the veins
- Refreshes + cools the legs
- Nourishes + moisturizes
- Medicinal plant extracts from horse chestnuts, sage, comfrey and witch hazel
- Care complex with evening primrose and apricot kernel oil

The tonic combines proven venous agents and sustainable care substances to create a unique and well-tolerated composition. In combination with support and compression stockings, blood circulation and microcirculation can be improved. Uncomplicated use! Can also be sprayed directly onto the stockings.

Ingredients: Aqua, Alcohol Denat, Glycerin, PEG-40 Hydrogenated Castor Oil, Sodium Lactate, Menthol, Lactic Acid, Hamamelis, Virginiana Leaf Extract, Litsea Cubeba Fruit Oil, Oenothera Biennis Oil, Prunus Armeniaca Kernel Oil, Symphytum Officinale Root Extract, Aesculus Hippocastanum Seed Extract, Camphor, Juniperus Oxycedrus Wood Oil, Linalool, Geraniol

COMPRESSANA VLIES Practical care tissues for skin and band

We make compression treatment easy and comfortable.

Compressana VLIES is used to care for areas of skin under compression treatment – even in the area of the band – and to clean the silicone coating if necessary. The fluid contained in the care tissue protects sensitive skin from irritation, soothes itching, moisturises and removes flaky skin. Skin irritations are prevented. It also ensures a better hold and fit of support stockings and prevents slipping. **Ingredients:** Aqua, Caprylic/Capric Triglyceride, Pentylene Glycol, Laureth-9, PEG-40 Hydrogenated Castor Oil



COMPRESSANA SKIN PROTECT Intensive care and repair creme

Skin regeneration and acute care of stressed, sensitive areas of the skin, even in cases of venous disorders. The antimicrobial effect of silver reduces the risk of infection. The repair complex can counter inflammation, irritation, etc. and promotes skin renewal.

For regular use, especially after wearing support and compression stockings, recommended especially on irritated skin areas. Also suitable for neurodermatitis and diabetics.

Ingredients: Aqua, caprylic/ capric Triglyceride, Panthenol, Glycerin, Pentylene Glycol, Polyglyceryl-3 Methylglucose Distearate, Zinc Oxide, Glyceryl Stearate, Stearyl Alcohol, Tocopherol, Silver, Citric Acid



Product	Art.No.	RRP
TONIC Spray 200 ml	0953	15,-€
TONIC Spray 50 ml	0952	6,50 €

Product	Art.No.	RRP
VLIES Care tissues in sachets, set of 10	0954	6,70 €

Product	Art.No.	RRP
SKIN PROTECT Tube 75 ml	0757	12,-€
SKIN PROTECT Mini-tube 10 ml	0756	6,-€

COMPRESSANA VEIN BALM

The vein balm protects and nourishes the skin and gives a pleasant and fresh feeling before and after wearing compression or support stockings.

Absorbs quickly and without residue. Selected natural ingredients (horse chestnut, arnica fluid extract and essential oils such as eucalyptus, camphor and rosemary) provide additional help and skin care for vein problems.

Ingredients: Aqua, Paraffinum Liquidum, Isopropyl Alcohol, Caprylic/Capric Triglyceride, Glycerin, Urea, Cetearyl Alcohol, Ammonium Acryloyldimenthyltaurate/Vp Copolymer, Trilaureth-4 Phosphate, Polyglyceryl-2-Sesguiisostearate, Eucalyptus Globulus Oil. Dimethicone. Camphor. Arnica Montana Extract. Aesculus Hippocastanum Seed Extract. Phenoxyethanol. Ethylhexylolycerin.



Product	Art.No.	RRP		
VEIN BALM Dash 125 ml	0104	11,50 €		
VEIN BALM Mini-tube 50 ml	0113	5,-€		

COMPRESSANA DERMA LOTION

Highly intensive emollient and skin protection for varicose veins, lipoedema, lymphatic oedema and skin diseases.

The lotion provides intensive and long lasting moisturisation. Protects the skin, keeps it supple

and resistant and prevents irritations and itching. For regular use after wearing compression, support and functional stockings.

Ingredients: Aqua, Isopropyl Palmitate, Caprylic/Capric Triglyceride, Butylene Glycol, Urea, Diisostearovl Polyglyceryl-3 Dimer Dilinoleate. Dimethicone and Dimethiconol Magnesium Sulfate, Hydrogenated Castor Oil, Cera Alba, Phenylpropanol, Tocopherol, Bisabolol, Backhousia Citriodora Leaf Extract Citral

> Bottle made of RECYCLED plastic

Product	Art.No.	RRP
DERMA LOTION Bottle 500 ml	0950	24,-€
DERMA LOTION Bottle 125 ml	0925	17,-€

COMPRESSANA HYDRO I OTION

For improved tolerance and better grip of stockings. Also suitable for application in the area of the grip top and silicone adhesive edge. before and after wearing compresion or support stockings.

Acute care against itching and for protection against irritations. Calms and soothes sensitive skin from mechanical irritations and rubbing. Enhances the skin regeneration and provides long lasting moisturisation.

Ingredients: Aqua, Caprylic/Capric Triglyceride, Glycerin, Pentylene Glycol, Polyalyceryl-3 Methylalucose Distearate, Laureth-9.Stearvl Alcohol. Glyceryl Stearate, Zinc Oxide, Tocopherol, Citric Acid, Xanthan Gum



Product	Art.No.	RRP
HYDRO LOTION Bottle 200 ml	0920	19,-€
HYDRO LOTION Airless-Spray 50 ml	0951	11,-€

Bottle

ACCESSORIES

Shoehorn long

This shoehorn makes it easier you to put on your shoes.

- Standing, with restricted mobility.
- with hook and practical hanging hole


Pocket nail file

Nail file in practical pocket style. Always ready to use and quickly to hand. Also ideal for on the go. To protect your stockings and for strong smooth nails, we recommend weekly nail care.

Double foot scraper

Double-sided callus file for gentle foot care that is kind to the skin. Pre-treat the callus with the coarse rubbing surface and smooth it with the fine side. Then apply SKIN PROTECT or a rich foot cream/callus balm. We recommend regular use to protect your stockings and for healthy, beautiful feet.







COMPRESSANA

Special bag – Tip for travelling and storage

Sharp edged or pointed objects can damage the knitting of your valuable stockings and tights In this special bag they can be kept protected and transported Helps you keep order at home and on journeys.



Product	Art.No.	RRP
Special bag for support and compression hosiery, 23 x 28 cm	0103	24,-€

Tip for the perfect care of the stockings

The special detergent COMPRESSANA CLEAN is suitable for elastic and high-quality textiles. Also for wool, silk and synthetic high-tech garments.

Clothes pegs – made of untreated beech wood



Product	Art.No.	RRP
CLEAN Special detergent Tube 250 ml	0121	4,90 €
CLEAN Special detergent Travel tube 50 ml,	0085	2,25 €
Clothes pegs made of beech wood, 6 pieces., 90 mm	0708	1,95 €



Laundry bag

Protects your stockings and tights during machine wash Quick and easy to fill by means of zipper.



Silver laundry bag

With pure silver fabric (0% nano particles) to combat unpleasant smelling laundry, even at low temperatures.

Product	Art.No.	RRP
Laundry bag with zipper, 100% Polyester	1007	4,50 €
Silver Laundry bag ith pure silver fabric, 50% PES, 50% PA-silver coated	0891	26,95 €

INSHOE-SOCKS



COMPRESSANA STEPS The tip for feet and stockings

For even better wearing comfort, flat elastic band at the end, not constricting, especially elastic and soft, air-permeable, not bulky, more comfortable than barefoot in the shoe.

Protect valuable stockings from damage in the foot area and create a healthy and pleasant climate in the shoe.



COMPRESSANA INSHOE SOCKS for protection and warm feet

Fashionable short shaft socks. Slipped over compression or support stockings/tights, the inshoesocks provide for warm feet during the cold season and protect the meshes against damage caused by footwear or little stones

Product	Colour	Size	Material	Art.No.	RRP
Inshoe-Socks with wool	anthrazit	35-38	80% wool, 20% PA	1003	8,-€
Inshoe-Socks with cotton	black, natur	39–42 43–45	98% cotton, 2% PA	1004	7,50€
STEPS with cotton	silk	35-38, 39-42	60% cotton, 40% PA	1005	4,-€



COMPRESSANA STANDS FOR HIGH QUALITY

Individuality:

In our range, you will always find the right sock in all different styles, colours and equipment or band variants.

Safety:

Our tested sizing system with six particularly tolerant standard sizes in three different lengths or exact custom-made sizes ensures optimal relief through a precisely defined pressure curve.

Guaranteed value:

Fine yarns and masterful manufacturing make for genuine quality stockings with pleasant wearing comfort and a beautiful stitch pattern.

Your quality and fit guarantee:

We will replace our compression products free of charge for 6 months from the date of purchase. The prerequisite is proper handling according to our instructions for use and care and the correct recording of body measurements or the correct determination of size.

Goodwill:

Despite proper use, damage to fine meshes due to unfortunate circum-stances, carelessness or external influences - who hasn't experienced this? Instead of "comprehensive insurance", we offer a goodwill arrangement for the "worst case scenario": replacement with new goods at a goodwill price within 6 months.

Sustainability:

COMPRESSANA compression stockings are harmless to human ecology and health and have been awarded the Oeko-Tex Standard 100 quality seal. Compressana does not use additives such as silicone or other softeners in the knitted fabric. All Compressana compression stockings are manufactured according to the highest quality standards, in an environmentally friendly and resource-saving manner in one of the most modern knitting factories in Germany.







Download the COMPRESSANA App on your tablet or smartphone. Use the new COLOURS 2021 consulting tool for easy colour selection of accessories and compression stockings. In addition, you will find selected leaflets and flyers and constantly supplemented and expanded information on our products.



App in German language availabe at Google Play Store and App Store.



COMPRESSANA GmbH Böhmerwaldstraße 3 · 93073 Neutraubling · Germany Tel.: 09401 9226-0 · Fax: 09401 9226-20 info@compressana.de · www.compressana.de